



Tosoh AIA Test Menu

Unit Dose Test Cup Reagent System

- Dry reagent format
- 90 day calibration stability for most assays
- Single, unitized cups - No pre-mixing or measuring
- Interchangeable results between Tosoh AIA Systems
- Bar-coded for identification/inventory management
- 10 minute incubation time for most assays



One Reagent for all Tosoh AIA Systems

The AIA-PACK Test Cup format works with every Tosoh AIA Automated Immunoassay System. From the compact AIA-360 to the full-featured AIA-600 II to the high-throughput AIA-2000, Tosoh's standardized Test Cup format ensures consistent performance every time.

Reagents are interchangeable, so laboratories with more than one Tosoh system can use the same reagents. Transition between systems is seamless, ensuring consistent results and efficient, economical operation for the laboratory.



AIA-360



AIA-600 II



AIA-2000

TOSOH BIOSCIENCE

Tosoh AIA Test Menu

10 minute incubation unless otherwise noted.

Tumor Markers

AFP
CEA
PSA
CA 125
CA 27.29
CA 19-9

Diabetes Assays

C-Peptide
HbA1c
Insulin

Reproductive

β -HCG
Estradiol
FSH
LH II
Prolactin
Progesterone
Testosterone

Metabolic

Cortisol
HGH

Thyroid

TSH
TT3
T4
FT3
FT4
T-Uptake
TSH 3rd Gen ♦
TPO Ab ♦
Tg Ab ♦

Cardiac Markers

CKMB
Troponin I 2nd Gen
Myoglobin

Anemia


Ferritin
Vitamin B12 ♦
Folate ♦
RBC Folate ♦

Kidney Marker Assays

β -2 Microglobulin
Cystatin C [○]
iPTH

Additional Assays

IgE II
PAP




Reagent Tray - 20 Test Cups/tray

[○] Pending CLIA categorization on AIA-360, AIA-600 II and AIA-2000

♦ 40 minute incubation. Available on the AIA-2000, AIA-600 II and AIA-1800 only

Unit Dose Test Cup Reagents

- Single, unitized cups use a dry reagent format that requires no pre-mixing, no pre-measuring, no on-board refrigeration and no waste.
- Dry reagent format ensures calibration stability of up to 90 days for most assays.
- Reagents have a one year shelf life from date of manufacture.
- The AIA test cups and trays are labeled with code and lot number for automated scheduling and inventory management on the AIA-2000 and AIA-1800.
- ST AIA-PACK reagents provide a faster turnaround time due to a 10 minute incubation protocol.
- Extensive test menu provides a wide range of immunoassay tests including thyroid, cardiac markers, tumor markers, reproductive hormones, anemia markers, metabolic markers and more.



Test Cup Diagram



Tosoh Bioscience, Inc.
6000 Shoreline Court, Suite 101
South San Francisco, CA 94080
Tel: (800) 248-6764
Fax: (650) 615-0415
www.tosohbioscience.us