



## miniVIDAS™

### The Compact Automated Immunoanalyzer

- Easy to use
- Cost-effective
- Reliable



BIOMÉRIEUX

# miniVIDAS™

## The Compact Automated Immunoanalyzer

### Easy to Use

- Always ready with no startup and minimal maintenance
- Easily operated by all personnel
- Ready-to-use reagents – just load and go

### Cost Effective

- Single test reaction strip
- Calibration curve provided for each test
- Each box contains everything needed to run 30 or 60 assays\*
- Batch or multiple assay testing

### Reliable

- Excellent CVs
- Most results in less than 30 minutes
- Reliable system with MTBF\*\* >1,000 days

\* To comply with CLIA guidelines, an additional control may need to be ordered.

\*\*Mean Time Between Failure



### Assays

Assay Name	Product Number	Tests Per Kit	Time (Minutes)	Reportable Range	CPT Code
<b>TUMOR MARKERS</b>					
TPSA	30428-01	60	60	0.07-100 ng/ml	84153
CEA-S	30453-01	60	60	0.5-200 ng/ml	82378
AFP	30413-01	60	30	0.6-400 ng/ml	82105
CA-125 II	30426-01	30	60	4-600 U/ml	86304
CA-15-3	30429-01	30	60	2-400 U/ml	86300
CA 19-9**	30427-01	30	60		86301
FPSA	30457	30	60	0.05-10 ng/ml	84153
<b>THYROID HORMONES</b>					
TSH	30400-01	60	40	0.05-60 µIU/ml	84443
FT4	30401-01	60	40	1-70 pmol/l	84439
T4	30404-01	60	40	6-320 nmol/l	84436
T3	30403-01	60	40	0.4-9 nmol/l	84480
<b>REPRODUCTIVE HORMONES</b>					
HCG	30405-01	60	30	2-1500 mIU/ml	84702
LH	30406-01	60	40	0.1-100 mIU/ml	83002
FSH	30407-01	60	40	0.1-110 mIU/ml	83001
Testosterone	30418	30	60	0.1-13 ng/ml	84403
Progesterone	30409-01	60	45	0.25-80 ng/ml	84414
Prolactin	30410-01	60	40	0.5-200 ng/ml	84146
Estradiol II	30431-01	60	60	9-3000 pg/ml	82670
<b>INFECTIOUS DISEASE</b>					
Measles IgG	30219	60	35	Qualitative	86765
Mumps IgG	30218	60	35	Qualitative	86735
Rubella IgG	30226	60	40	0-250 IU/ml	86762
CMV IgG	30204-01	60	40	0-400 AU/ml	86644
CMV IgM	30205-01	30	60	Qualitative	86645
Toxo IgG II	30210-01	60	40	0-300 IU/ml	86777
Toxo IgM	30202-01	60	40	Qualitative	86778
Toxo Competition	30211-01	60	40	Qualitative	86777
Varicella Zoster IgG	30217	60	35	Qualitative	86787
Chlamydia*	30101-01	60	60	Qualitative	87320
Chlamydia Blocking*	30194-01	30	75	Qualitative	87420
<i>C. difficile</i> Toxin A&B	30118-01	60	75	Qualitative	87324
Rotavirus*	30107	60	30	Qualitative	86759
<b>THERAPEUTIC DRUG MONITORING</b>					
Digoxin	30603-01	60	20	0.2-5.0 ng/ml	80162
<b>SEROLOGY</b>					
<i>Helicobacter pylori</i> IgG*	30192	30	35	Qualitative	86677
Lyme IgG/IgM	30298-01	60	35	Qualitative	86618
<b>COAGULATION</b>					
D-Dimer Exclusion®	30442	60	35	45-10,000 ng/ml	85379
<b>CARDIAC ASSAYS</b>					
Troponin I Ultra	30448-01	60	20	0.01-30 ng/ml	84484
CKMB	30421-01	30	35	0.8-300 ng/ml	82553
NT-proBNP	30449-01	60	20	20-25,000 pg/ml	83880
<b>OTHER ASSAYS</b>					
Ferritin	30411-01	60	30	1.5-1200 ng/ml	82728
Total IgE	30419-01	60	30	0.5-1000 kIU/l	82785
Procalcitonin PCT	30450-01	60	20	0.05-200 ng/ml	0194T

\*Highly complex

\*\*Pending FDA clearance/approval, not available for sale or distribution in the U.S.

bioMérieux, Inc.  
100 Rodolphe Street  
Durham, NC 27712  
U.S.A.

Tel: (800) 682 2666  
Fax: (800) 968 9494

bioMérieux Canada, Inc.  
7815 Henri-Bourassa Boulevard West  
St. Laurent (Québec)  
Canada H4S 1P7

Tel: (1) 514 336 7321  
Fax: (1) 514 807 0015

[www.biomerieux-usa.com](http://www.biomerieux-usa.com)

